

Assays-in-a-Box™ by Preferred Cell Systems™

HALO QC[®]

A Stem Cell Quality Control Assay for Bone Marrow, Umbilical Cord Blood and Peripheral Blood Cellular Therapy Processing

Technical Manual

(Version 10-23)

This manual should be read in its entirety prior to using
this product

For *In Vitro* Research Use Only.

Not for clinical diagnostic use.

No part of this instruction manual may be copied, duplicated or used
without the express consent of Preferred Cell Systems™

Preferred Cell Systems™

TABLE OF CONTENTS

1. Limitations of the Assay and Precautions	3
2. Introduction	4
3. Use and Availability	4
4. The Concept and Principle of ATP Bioluminescence Assays	5
5. QuickGuide to HALO® QC	7
6. Kit Contents and Storage Conditions	8
7. Equipment, Supplies and Reagents Required, but not Provided	9
8. The HALO® QC Protocol	9
Step 1 - Cell Preparation	9
Step 2 - HALO® QC Cell Culture Procedure	11
Step 3 - Bioluminescence Measurement	12
9. Recommendations and Tips Prior to Using HALO® QC	13
10. Recommendations and Tips Prior to Measuring Bioluminescence	14
11. Luminescence Plate Reader Setup and Conversion of RLU Values to ATP Values using the ATP Standard Curve	16
12. Results	16
13. HALO® QC Assay Measurement Assurance and Validation Parameters	17
14. Troubleshooting	18
Other Hematopoietic Cellular Therapy Assays Products	21
Calibration and Assay Standardization - Protocol 1	23
Calibration and Assay Standardization - Protocol 2	24
KM016.002	

1. LIMITATIONS OF THE ASSAY AND PRECAUTIONS

1. HALO® QC is not approved by either the U.S. Food and Drug Administration (FDA) or the European Medicines Agency (EMA)
2. HALO® QC is for research use only and has not been approved for clinical diagnostic use.
3. Reagents and supplies in this kit are STERILE. Perform all procedures under sterile conditions, except where indicated.
4. This kit should not be used beyond the expiration date on the kit label.
5. Do not mix or substitute reagents or other kit contents from other kit lots or sources.
6. Always use professionally calibrated and, preferably, electronic pipettes for all dispensing procedures. Small discrepancies in pipetting can lead to large pipetting errors. Although electronic pipettes self-calibrate themselves, they still need to be professionally calibrated on a regular basis.
7. Good laboratory practices and universal protective precautions should be undertaken at all times when handling the kit components as well as human cells and tissues. Material safety data sheets (MSDS) are included in each literature packet.

2. Introduction

Present quality assurance for umbilical cord blood, bone marrow or mobilized peripheral blood cellular therapy products is based on four parameters, namely (a) total nucleated cell (TNC) count, (b) dye exclusion viability, (c) CD34 content by flow cytometry and (d), the colony-forming unit (CFU) assay. None of these tests can provide information on the quality of the stem cells in the original or processed product that would indicate loss or damage to these important cells required for engraftment and reconstitution.

Stem cell quality is dependent upon a number of factors, including how the cells are obtained, how the tissues are processed to obtain the stem cell product and whether the cells are cryopreserved. At each stage of the process to manufacture and produce a stem cell product, a sample should be tested to (a) ensure viability, stem cell number and yield, (b) ensure manufacturing procedure consistency and (c), ensure stem cell proliferation ability. Only by measuring stem cell quality using a standardized and validated assay at each procedure step is it possible to ensure a high quality product.

HALO® QC (Quality Control) was designed specifically for this purpose. HALO® QC is used to detect either (a) the primitive hematopoietic stem cell population, SC-GEMM 1 (equivalent to colony-forming CFC-GEMM, granulocyte, erythroid, macrophage, megakaryocyte), or (b) both the SC-GEMM 1 and an even more primitive lympho-hematopoietic stem cell population designated SC-HPP 2 (which does not have a CFU equivalent).

Stem cell proliferation ability is measured by virtue of the fact that when stem cells are stimulated with a growth factor cocktail for each of these stem cell populations, the intracellular ATP (iATP) concentration increases proportionately. Indeed, the iATP concentration and, therefore the proliferation status, correlates directly with the number of stem cells plated. Therefore, to obtain a measure of stem cell quality, the iATP concentration of the cell population is detected at a single cell dose.

HALO® QC incorporates a fully standardized and validated ATP bioluminescence readout to measure stem cell proliferation ability. HALO® QC is, therefore, an instrument-based, non-subjective measurement of iATP using a luciferin/luciferase reaction. It is the most sensitive and accurate, non-radioactive readout available. HALO® QC produces reliable and reproducible results, not only within a single laboratory, but between different laboratories.

Like all HALO® assays, HALO® QC also incorporates Suspension Expansion Culture™ (SEC) Technology. SEC™ technology is a methylcellulose-free culture system that is fast and easy to use, exhibits greater sensitivity, accuracy and precision than any CFU assay and requires significantly shorter culture times than the CFU assay. The investigator can be assured that when using HALO® QC in the stem cell processing laboratory, they will obtain reliable and trustworthy results that can be compared over time.

3. Use and Availability

HALO® QC is intended for use in hematopoietic stem cell processing laboratories that manufacture and produce hematopoietic stem cell products for transplantation purposes. HALO® QC should be used to determine stem cell quality (proliferation ability) at each stage of the manufacturing process.

HALO® QC is used to determine the quality of stem cells derived from human:

- Umbilical cord blood
- Bone marrow
- Peripheral blood (normal or mobilized)

HALO® QC is available for measuring the stem cell quality of:

- *In vitro*, primitive hematopoietic stem cells, designated SC-GEMM 1 (equivalent to CFC-GEMM or Colony-Forming Cell - Granulocyte, Erythroid, Macrophage, Megakaryocyte).

Preferred Cell Systems™

- Both the SC-GEMM 1 and the in vitro, primitive lympho-hematopoietic stem cell population, designated SC-HPP 2 (Stem Cell- High Proliferative Potential 2).

Although HALO® QC can be used to detect stem cells in a total nucleated cell (TNC) fraction, it is NOT recommended since (a) a TNC fraction contains cell impurities that dilute the stem cells, (b) this results in an underestimation of stem cell proliferation ability, and (c) interference by high red blood cell concentrations with the ATP readout results in a false positive result.

The cell preparation used for HALO® QC should be a:

- Mononuclear cell (MNC) fraction.
- Purified stem cell fraction, e.g. CD34+.

Suspension Expansion Culture™ (SEC™) Technology

HALO® QC incorporates Suspension Expansion Culture (SEC) Technology. No methylcellulose is used. This has the following advantages over traditional CFU assay methodology:

- All reagents can be dispensed using normal pipettes.
- Allows cell-cell interaction.
- Produces greater assay sensitivity.
- Shorter cell incubation times; cell proliferation is measured on the exponential part of the growth curve.
- Coefficients of variation $\leq 15\%$.

HALO® QC can be obtained with serum-free or low serum SEC™ Master Mix formulations.

HALO® QC Assay Availability

Catalog No.	Cell Populations(*)	Serum Formulation	No. of Samples	No. of Plates
K2-1QC-1	SC-GEMM 1	Low serum	16-24 samples	1
K2SF-1QC-1	SC-GEMM 1	Serum-free	16-24 samples	1
K2-2QC-2	SC-HPP 2 + SC-GEMM 1	Low serum	16-32 samples	2
K2SF-2QC-2	SC-HPP 2 + SC-GEMM 1	Serum-free	16-32 samples	2

(*) To distinguish between colony-forming cells (CFC) and cells that are grown in suspension cultures (e.g. HALO®), stem cells are designated with the prefix, SC.


4. The Concept of ATP Bioluminescence Assays

HALO® QC is an ATP bioluminescence assay. The fundamental concept underlying the assay is the measurement of the cell's chemical energy in the form of intracellular adenosine triphosphate (iATP). If a cell is producing iATP, it is demonstrating cellular and mitochondrial integrity and is therefore viable. When hematopoietic cells are stimulated to proliferate, in culture with growth factors and/or cytokines, the iATP concentration increases several fold. The iATP concentration produced is directly dependent on:

- The proliferation potential (or primitiveness) of the cell population being detected.
- The concentration of the growth factor(s) and/or cytokine(s) used to stimulate the cells.
- The plated cell concentration.

Preferred Cell Systems™

Hematopoietic stem cells are incubated in the HALO® QC Master Mix provided with this kit for a specific period of time. When the culture period has elapsed, a single-step addition of an ATP Enumeration Reagent is dispensed into each culture well and the contents mixed. The plate is incubated at room temperature in the dark for 10 minutes. During this time, the cells are lysed and the released iATP acts as a limiting substrate of a luciferin/luciferase reaction to produce bioluminescence in the form of light according to the following equation:



The bioluminescence emitted is detected and measured in a luminescence plate reader as relative luminescence units (RLU). The assay can be calibrated and standardized, and controls and standards are included for this purpose. Performing an ATP standard curve and controls has the following advantages:

1. The controls calibrate the instrument and also ensure that the reagents are working correctly.
2. The ATP standard curve also ensures that the reagents are working correctly.
3. The ATP standard curve allows the luminometer output in Relative Luminescence Units (RLU) to be converted to standardized ATP concentrations (μM).
4. Performing the ATP standard curve allows results to be compared over time.
5. The results obtained from controls and standard curve should be compared with those provided in Section 12. These are the measurement assurance parameters that allow the investigator to ensure that the assay is working correctly prior to measuring samples. When the values from the controls and ATP standard curve are within the ranges provided in Section 12, the investigator can consider the results trustworthy.

The ATP standard curve and controls need only be measured once on the day samples are to be processed. Do not use previous results from an ATP standard curve and controls performed on a different day. This will cause erroneous results.


The ATP standard curve is used to convert sample RLU results into ATP concentrations by interpolation. This procedure can often be performed automatically by the instrument software. If the software does not allow this, it will be necessary to use third-party software to perform this operation.

NOTES

5. QuickGuide to HALO® QC (Figure 1)

QuickGuide to HALO® QC

STEP 1
 If possible, fractionate cells to a mononuclear cell (MNC) fraction as this will provide the best results.
 Perform cell count and viability and adjust to working cell concentration.
 A minimum of 6 replicates wells/sample is recommended.



QG01 804

6. Kit Contents and Storage Conditions

HALO® QC kits contain reagents that have been frozen and stored at -80°C prior to shipment. The kit is shipped either with dry ice or blue ice. The following components are included:

Item	Component	Storage
1	HALO® QC Master Mix to detect primitive hematopoietic stem cells (SC-GEMM 1).	-20°C until used
1a	HALO® Master Mix to detect primitive lympho-hematopoietic stem cells (SC-HPP 2) and SC-GEMM 1.	-20°C until used
2	Medium (IMDM) for dilution of the ATP standard.	-20°C until used
3	ATP standard.	-20°C until used
4	ATP extra high, high and low controls.	-20°C until used
5	ATP Enumeration Reagent (ATP-ER)*	-20°C in the dark until used
6	Adhesive Plate Covering: a sterile foil to protect and keep unused wells sterile.	Can be kept with other kit components
7	Sterile, solid white 96-well plates for stem cell culture	Can be kept with other kit components
8	Non-sterile, solid white 96-well plates for ATP standard curve determination.	Can be kept with other kit components
	Technical manual can be downloaded from the HALO® QC page on the Preferred Cell Systems™ website	

Exact volumes of the kit reagents and supplies are provided on a separate sheet included with this assay kit.

*The ATP-ER should not be thawed until needed and can be refrozen 11 times without significant loss of sensitivity. It can be kept at 2-8°C for 48h once thawed and is stable for 20 weeks when stored at -20°C. This reagent is light sensitive. Keep in the dark. The ATP-ER must not be used past the expiration date.

IMPORTANT

All kit components are quality controlled and optimized so that they work together. Please do not replace kit components with those of a different product. This will invalidate the warranty provided by Preferred Cell Systems™.

This kit contains a reagent for measuring luminescence (ATP-ER) that decays with time. Preferred Cell Systems™ recommends that this kit be used before the expiry date of this reagent. Preferred Cell Systems™ does not take responsibility for the quality of reagents beyond their expiry date. If the kit cannot be used prior to the expiry date of this reagent, fresh reagent can be obtained from Preferred Cell Systems™.

Good laboratory practices and universal protective precautions should be undertaken at all times when handling the kit components as well as cells and tissues. Material safety data sheets (MSDS) are included in each literature packet.

7. Equipment, Supplies and Reagents Required, but not Provided

Equipment and Supplies

1. Laminar Flow Biohood
2. Plate luminometer (e.g. Berthold LB962 CentroLIA/pc; Molecular Devices, SpectraMaxL)
3. Sterile plastic tubes (5ml, 10ml, 50ml)
4. Single channel pipettes, preferably electronic (e.g. Rainin EDP pipettes for variable volumes between 1µl and 1000µl)
5. 8 or 12-channel pipette, preferably electronic (e.g. Rainin EDP pipettes for fixed or variable volumes between 10µl and 100µl)
6. Reservoir for 8- or 12 channel pipette
7. Sterile pipette tips
8. Vortex mixer
9. Tissue culture incubator, humidified at 37°C with 5% CO₂ (minimum requirement) and 5% O₂ (preferable)
10. 1.5ml plastic vials (5 for each ATP dose response)
11. Hemocytometer or electronic cell counter to determine cell concentration
12. Flow cytometer or hemocytometer for determining viability.

Reagents

1. HemoGro™ Hematopoietic Growth Medium for cell suspensions and dilutions (Preferred Cell Systems™)
2. Iscove's Modified Dulbecco's Medium (IMDM)
3. Density-gradient medium (e.g. Ficoll-Paque, Lymphoprep).
4. 7-AAD, propidium iodide or trypan blue for viability assay.
5. LIVEGlo™ metabolic viability assay (Preferred Cell Systems™)

8. The HALO® QC Protocol

**PLEASE READ THE FOLLOWING PROTOCOL CAREFULLY.
SEE SECTION 9 BEFORE PERFORMING THE PROTOCOL.**

Performing HALO® QC is a 3-step process.

Step 1 – Cell preparation.

Step 2 – HALO® QC cell culture master mix procedure, plating and incubation in the 96-well plate.

Step 3 – Luminescence measurement. An ATP dose response is performed prior to sample luminescence measurements with conversion of RLU to µM ATP.

Step 1 and Step 2 must be performed in a laminar flow biohazard hood

STEP 1 – Cell Preparation

HALO® QC can be performed using tissues with the following purity:

1. Total nucleated cell (TNC) fraction usually produced by red blood cell reduction. The TNC fraction is often used for human bone marrow and umbilical cord blood. The concentration of red blood cells in this preparation may be 30% or higher. Although the TNC fraction can be used, it is NOT recommended due to (a) dilution of primitive stem and progenitor cells, (b) underestimation of primitive stem and progenitor cells, and (c) interference of high red blood cell concentrations with

Preferred Cell Systems™

the ATP readout resulting in an unsatisfactory high ATP readout. It is highly recommended to use a mononuclear cell fraction or higher for these assays.

2. Mononuclear cell (MNC) fraction is the cell preparation of choice. This fraction can be prepared by density gradient centrifugation using, for example, Ficoll-Paque or Lymphoprep.
3. Purified stem and progenitor cell populations prepared by cell sorting or magnetic bead separation.

Cell Viability, Cell Counting and Cell Culture Suspension Preparation

1. For dye exclusion viability methods, use trypan blue and a hemocytometer or automated method, flow cytometry using 7-AAD or another vital stain.

Note that dye exclusion viability methods detect membrane integrity. They do not detect cellular and mitochondrial integrity and therefore metabolic viability.

A viability of 85% or greater should be obtained when using dye exclusion viability methods only. It is recommended not to use cell suspensions with a viability of less than 85% since these cells will not be able to sustain proliferation ability. Use LIVEGlo™ (Preferred Cell Systems™) as a metabolic viability assay.

2. Determine the cell concentration using either a hemocytometer or electronic cell/particle counter. **NOTE:** Do not base the working concentration on the TNC count or the number of viable cells as this will give erroneous results.
3. Adjust the cell suspension concentration to that recommended in Table 1.
Note the working cell concentration per ml is 100 x the final cell concentration per well. If cells have been treated prior to cell culture, higher cell concentrations may be required.
4. Prepare the total volume of cell suspension required using HemoGro™. The volume of the adjusted cell suspension required will be 10% of the total volume of HALO® QC Culture Master Mix prepared.

TABLE 1
Recommended Cell Doses for Cell Types, Cell Preparations and Cell States for HALO® QC

Species	Cell Type	Cell Preparation	Cell State	Working Cell Concentration Required (100 x Final Cells/Well)	Final Cell Dose / Well
Human	Bone marrow	MNC	Fresh/Frozen	0.5-0.75 x 10 ⁶	5,000-7,500
	Peripheral blood	MNC	Fresh/Frozen	0.5-0.75 x 10 ⁶	5,000-7,500
	Mobilized peripheral blood	MNC	Fresh/Frozen	0.5-0.75 x 10 ⁶	5,000-7,500
	Umbilical cord blood	MNC	Fresh/Frozen	0.5-0.75 x 10 ⁶	5,000-7,500
	Umbilical cord blood	MNC	Frozen	0.5-0.75 x 10 ⁶	5,000-7,500
	Bone marrow	CD34 ⁺	Fresh	0.1-1 x 10 ⁵	100-1,000
	Mobilized peripheral blood*	CD34 ⁺	Fresh/Frozen	0.1-5 x 10 ⁵	100-5,000
	Umbilical cord blood	CD34 ⁺	Fresh/Frozen	0.1-5 x 10 ⁵	100-5,000

STEP 2. HALO® QC Cell Culture Procedure

- HALO® QC Master Mixes are complete and ready-to-use.
- Perform all procedures under a laminar flow, bio-hazard hood.
- Wear protective clothing, including gloves for all operations.

The HALO® QC Methodology

1. Remove the HALO® QC Master Mix that came with the assay kit from the freezer and thaw either at room temperature or in a beaker of cold water. Do not thaw in a 37°C water bath or incubator.
2. Label sufficient 5mL tubes for the number of samples to be tested.
3. Using a calibrated pipette, preferably an electronic pipette with repeat function, dispense a volume of the HALO® QC Master Mix minus 10%, that will allow 0.1mL to be dispensed into the required number of replicate wells. See Table 2.

TABLE 2

Number of Replicate Wells Required	Volume of HALO® QC Master Mix	Volume of Cells (10% of final volume)	Total Volume
6	0.63mL	0.07mL	0.7mL

4. Adjust the cell concentration to the working cell concentration shown in Table 1. The working cell concentration should be 100 fold higher than the final cell concentration in the well. Thus, if the final cell concentration is to be 5,000 cells/well, the working cell concentration should be 100 times 5,000 or 500,000 (5×10^5) cells/mL.
5. Dispense the required volume of cells into each tube containing the HALO® QC Master Mix.
6. Mix the contents gently on a vortex mixer. Do not cause bubbles.
7. Remove the sterile, individually wrapped, 96-well plate from the assay kit box.
8. Using a calibrated pipette, preferably an electronic pipette with repeat function, dispense 0.1mL into each of the required number of replicate wells of the 96-well plate.
TIP: If preparing 4 or 6 replicates, dispense either across the plate in rows or if 8 replicates are prepared, these should be dispensed across the plate in columns.
9. After replacing the lid, transfer the 96-well plate to a humidified container (see Section 9).
10. Transfer the humidified container to the incubator. The incubator should be fully humidified at 37°C and gassed with 5% CO₂ and, if possible, 5% O₂. This reduced oxygen tension improves plating efficiency by reducing oxygen toxicity caused by the producing of free radicals.
11. Incubate the cells for 5 days. For greater sensitivity, the incubation time can be extended to 7 days. The ATP concentration should increase 2-3 fold from day 5 to day 7. However, this may occur at the expense of slightly greater variation between the replicates.

STEP 3 – Bioluminescence Measurement

Please note the following important points:

- FOR ALL OF THE FOLLOWING STEPS, WEAR LABORATORY GLOVES. ATP is present on the skin and can cause erroneous results
- PLEASE DOWNLOAD THE DOCUMENT ON HOW TO SETUP THE PLATE LUMINOMETER. The instrument should be setup and prepared for use prior to any of the following steps being performed.
- Please refer to Section 10 for recommendations and tips prior to starting this part of the procedure. In particular, please refer to Section 10 for important information on mixing components.
- Remove the ATP Enumeration Reagent (ATP-ER) from the freezer and thaw at room temperature or in cold running water prior to analysis. Do not thaw the ATP-ER in a water bath or 37°C incubator.
- If the assay is to be calibrated and standardized, remove the ATP standard, controls and reagents from the freezer and thaw to room temperature or in cold running water prior to analysis.
- ATP standard curves performed on previous days or for previous experiments or studies must not be used since the ATP-ER intensity changes with time and lot number.
- Use the unwrapped, non-sterile, 96-well plate provided with the kit to perform the ATP standard dose response curve.

A. Instrument Calibration and Assay Standardization

An Instructional Tutorial on “How to Calibrate and Standardize Any ATP Bioluminescence Assay” is available on the Preferred Cell Systems website under the Resources tab.

It is highly recommended to calibrate the luminescence plate reader and standardize the assay prior to measuring any samples. This will allow a comparison with the expected measurement assurance values (**see Section 13**) that should be obtained prior to measuring samples. Use the non-sterile, 96-well white plate provided with the assay kit for this purpose.

HALO® QC Assay Kit includes the following to calibrate and standardize the ATP bioluminescence part of the assay to measure cell proliferation.

- IMDM medium: Used only for ATP standard serial dilution.
- ATP Standard at 10µM: Serially diluted to produce the ATP standard curve.
- Low ATP Calibration Control: Used for normal and extra high cell proliferation.
- High ATP Calibration Control: Used for normal cell proliferation.
- Extra High ATP Calibration Control: Used for extra high cell proliferation.

B. Deciding Which Calibration Controls to Use and ATP Standard Curve Range

PROTOCOL 1: For bone marrow and umbilical cord blood, use the low and high calibration controls and perform an ATP standard curve from 0.01µM to 1µM. **See Page 23.** ATP concentrations should be within this ATP standard curve range. If they are found to be higher, perform Protocol 2.

PROTOCOL 2: For most mobilized peripheral blood products and purified cells (e.g. CD34⁺ cells), it is recommended to use the low and extra high calibration controls and perform an ATP standard curve from 0.03µM to 3µM. **See Page 24.** It should be noted that not all patients mobilized in the same manner and some may not mobilize at all. In this case, use Protocol 1.

If calibration and assay standardization has been performed, it is extremely important that the ATP standard curve slope, and ATP controls lie within the ranges specified in **Section 13**.

It is important that the sample ATP values are within the limits of the ATP standard curve, otherwise the interpolation of Relative Luminescence Unit (RLU) values from the luminescence plate reader into ATP concentrations will not be accurate. In some cases, cell proliferation could be greater than 3µM ATP. If ATP values from the samples are greater than 3µM, it is recommended to dilute the sample with additional medium so that the values are within the ATP standard curve range. This may require removing an aliquot from the replicate wells, transferring the aliquot to a new wells and diluting each aliquot with additional medium. The replicate wells would then be reread.

C. Sample Measurement

The addition of ATP-ER is performed in the same manner as the ATP Standard Curve.

1. If possible, place the sample plate(s) in a humidified incubator set at 22-23°C gassed with 5% CO₂ for 30min. Otherwise, allow the plate to come to room temperature for 30 min.
2. If only part of the plate has been used, transfer the plate to a bio-safety hood and remove the lid under sterile conditions. Take a sterile adhesive plate coverfoil from the kit box, remove the backing and layer it over the top of the plate. Using a sharp knife or scalpel, cut away the foil that covers the wells to be processed. The unused, empty wells will now remain sterile for the next samples. (See Section 10, Adhesive Plate Covering Film).
3. Using a multichannel pipette (8- or 12-channel depending on the plate configuration), add 0.1mL of ATP-ER to each well of the first column (A1-H1) or row (A1-12). Mix the contents as described in Section 10.
4. Repeat this procedure for each column or row using new tips.
5. When ATP-ER has been added to all wells, replace the plastic cover and incubate for 10 min at room temperature in the dark to lyse the cells and stabilize the luminescence signal. Incubate the plate in the reader for the last 2 min to stabilize the plate.
6. Unused ATP-ER may be returned to the bottle and refrozen. See section 10 for ATP reagent storage conditions and stability.

D. Using a plate luminometer with automatic dispenser

The user may have a plate luminometer that allows reagents to be dispensed automatically directly into the well. Preferred Cell Systems™ does **not** recommend using the automatic dispensers, since the contents of the well are not mixed sufficiently using this method.

9. Recommendations and Tips Prior To Using HALO® QC.

(i) Cell Suspension

- a) The preferred cell suspension is a mononuclear cell suspension (MNC) or higher purity. **PLEASE NOTE:** Although it is normal for cell processing laboratories to prepare a red blood cell-/plasma-depleted total nucleated cell (TNC) fraction using instrumentation, using this fraction for determining stem cell quality will be severely limited and could result in erroneous conclusions about the quality of the sample and the unit of cells in general.
- b) Extraneous ATP, red blood cells (which have high concentrations of ATP) and hemoglobin interfere with the ATP analysis. The cell suspension must have a hematocrit of 10% or less.
- c) If cells have been treated prior to cell culture, higher cell concentrations than those shown in Table 1 may be required.

(ii) Number of Replicates Performed

Six (6) replicates/samples are used for HALO® QC. Please remember that using fewer replicates may save components in the short term, but may also cause inconclusive results. If outliers are encountered, which may have to be removed from the analysis, the consequence could be that extra experiments would be required resulting in extra time and costs.

Preferred Cell Systems™

(iii) Plate Configuration

For 6 replicate wells/sample dispense in rows across the plate.

(iv) 96-Well Plates Provided

The reagents have been optimized to work with the 96-well plate(s) provided in the HALO® QC kit. Please do not replace the plates included with the kit with those of another manufacturer. Cell growth and bioluminescence output can be seriously affected and the assay kit warranty will be void. Additional plates can be purchased from Preferred Cell Systems™ if required.

(v) Humidity Chamber

A humidity chamber is recommended due to the small sample volume. Even fully humidified incubators do not keep the humidity level high enough to keep the sample from evaporating. This usually results in so-called “edge effects”. This phenomenon is observed when ATP values in the outside wells are lower than those in the inside wells. A humidity chamber can be assembled using plastic lunch boxes or other plasticware available from a supermarket or discount stores. Holes must be made in the lid to enable adequate gas exchange. Disposable serological pipettes are cut to an appropriate length to fill the bottom of the container. Distilled/deionized water is poured into the container to just below the level of the pipettes. This allows for adequate water to keep the humidity high without the plates sitting in water. Please contact Preferred Cell Systems™ for further information about assembling and using humidity chambers.

(vi) Incubation Times

The incubation time may vary depending on cell type and species. Assay sensitivity might improve with longer incubation times, but usually at the expense of higher variability between wells. Once an optimal incubation time has been found, the same time period should be maintained for all future experiments so that results can be directly compared.

10. Recommendations and Tips Prior To Measuring Bioluminescence

- *Always wear laboratory (e.g. latex) gloves during this operation to avoid ATP contamination from skin.*
- *DO NOT wipe the pipette tip with tissue etc as this will wick the reagent from the tip and cause an erroneous ATP standard curve and false sample results.*
- *Always change pipette tips after each use.*
- *Each day bioluminescence is measured, a standard curve MUST be performed. The ATP-ER decays with time. A new ATP standard curve must be performed to ensure accurate conversion of the RLU values to ATP concentrations so that results can be compared.*
- *HALO® QC includes solid white plates for both cell culture and the ATP standard curve and controls. Do not use different plates for the assay. Doing so will result in inaccurate results and invalidation of the assay kit warranty. Extra plates can be purchased from Preferred Cell Systems™.*

Bioluminescence Assay Kit Components

- Prior to measuring bioluminescence, remove the ATP standard, 1 set of ATP controls and the ATP-Enumeration Reagent (ATP-ER) from the freezer and thaw at room temperature or at 22 - 23°C.
- Sufficient ATP standard, controls and ATP-ER are supplied to usually perform 2 standard curves and controls/assay kit. Additional ATP standards and controls can be obtained from Preferred Cell Systems™.
- If thawing more than one bottle of ATP-ER for analysis, mix the contents of the bottles together before dispensing into reagent reservoir.

Preferred Cell Systems™

- ATP-ER can be refrozen up to 11 cycles without significant loss of sensitivity. Thawed ATP-ER can be kept at 2-8°C, in the dark, for 48h or is stable at -20°C for 20 weeks.

Reconstitution of Lyophilized Monitoring Reagent (if included)

- Thaw the ATP Enumeration Reagent Buffer at room temperature, in cold running water, or at 2-8°C overnight.
- Do not use any form of heat to thaw this reagent.
- Allow the lyophilized ATP-ER substrate (brown glass bottle) to come to room temperature.
- Remove the closures from both bottles.
- Carefully pour the entire contents of the buffer bottle into the lyophilized ATP-ER substrate bottle. Swirl gently or invert slowly to mix. Do not shake.
- Allow the ATP-ER mix to reconstitute for 10 minutes at room temperature.
- Reconstituted ATP-ER is stable for 8 hours at room temperature, 48 hours at 2-8°C, or 20 weeks at -20°C.
- ATP-ER can be refrozen up to 11 cycles without significant loss of sensitivity.

Volumes of Luminescence Kit Components Required

- Each vial of ATP standard contains enough volume to perform one or two ATP standard dose responses.
- The amount of ATP-ER added to each well is 0.10mL. Therefore:
Total amount of ATP-ER (µl) required = 0.1mL x (number of wells used + 24 (background, ATP dose response wells and ATP controls)).

ATP Standard Curve

Depending on the size of the kit purchased, non-sterile, 96-well plates have been included to perform an ATP standard curve prior to processing the sample cultures. Performing an ATP standard curve and controls on each day samples are processed is an essential part of the assay because it has 4 functions:

- It tests whether the instrument is working properly and calibrates it.
- It ensures that the reagents are working correctly.
- It calibrates and standardizes the assay and allows the assay system to be validated, if required.
- It allows the output of the plate luminometer, in relative luminescence units (RLU), to be converted to ATP concentrations, thereby standardizing the procedure so that intra- and inter-laboratory experiments can be compared.

Adhesive Plate Covering Film

To help keep the plate(s) sterile, adhesive, air permeable, sterile films are provided so that the part of the plate that is not being used can be covered and kept sterile until required. If using the adhesive film provided, the plate cover should be removed in a laminar air-flow hood and replaced with the film to ensure sterility.

Mixing the Contents of 96-well Plate


Please refer to the Instructional Tutorial on “How to Calibrate and Standardize Any ATP Bioluminescence Assay” is available on the Preferred Cell Systems website under the Resources tab.

Mixing the contents of the wells after adding ATP-ER is one of the most important procedures of the assay. It is recommended that the addition of ATP-ER is performed using a multi-channel pipette to achieve consistency and reduce variability. Addition of the reagent and mixing should be performed in the following manner:

1. Take up the required amount of reagent and add it to the well without inserting the tip into the well contents.
2. Starting from the center of the well, aspirate and dispense the contents twice without removing the pipette tip from the contents of the well.

3. Move the pipette tip to one corner of the well and aspirate and dispense the contents twice without removing the tip from the contents of the well.
4. Repeat this operation as shown in Figure 4 for each corner of the well.
5. Try not to cause excessive bubbles in the culture and **DO NOT** over mix since this can result in drastically reduced luminescence values.
6. This procedure effectively and optimally mixes the contents well.

Figure 2. Positions of pipette tip for mixing the well contents



11. Luminescence Plate Reader Setup and Conversion of RLU Values to ATP Values Using the ATP Standard Curve

It is very important that the luminescence or multimode plate reader is setup correctly, otherwise false results could occur. Preferred Cell Systems™ has provided a separate document to help the investigator setup their instrument and perform the calculations in order to convert Relative Luminescence Units (RLU) into ATP concentrations using the ATP standard curve. It is strongly recommended that the investigator consult this document prior to performing any ATP bioluminescence assay. **This document can be downloaded with this manual from the Preferred Cell Systems™ website.**

12. Results

HALO® QC provides an instrument-based, non-subjective, quantitative readout of metabolic viability and the total proliferation ability of the cells being tested. HALO® QC is an alternative to the colony-forming unit (CFU) assay for hematopoietic cellular therapy processing laboratories. Since HALO® QC measures cell proliferation and proliferation occurs prior to differentiation, HALO® QC can not only be completed more rapidly, but actually predicts results obtained with the CFU differentiation assay. This is because there is a direct correlation between HALO® QC and the CFU assay.

HALO® QC can be used on fresh or frozen cells. If testing cells prior to cryopreservation, it should be emphasized that results produced using HALO® QC or any other QC assay from Preferred Cell Systems™, will not reflect the actual viability or proliferation ability after thawing; proliferation ability will be 2-3 fold lower after thawing a cryopreserved sample than a fresh sample.

It is strongly recommended that the user perform the calibration and ATP standardization procedure described in this manual and available as an Instructional Video on the Preferred Cell Systems website under Resources. The ATP controls calibrate the luminescence plate reader. The ATP standard curve allows non-standardized RLU values to be converted into standardized ATP concentrations (μM). The results should be compared with those provided in **Section 13** below. Providing the results are within the ranges specified in **Section 13**, it is then possible to continue with sample processing and measurement.

These measurement assurance parameters also indicate whether the sample being tested provides acceptable viability and proliferation activity, or whether it should be rejected for use. Please be aware that acceptance criteria for the sample may not necessarily apply to the unit of cells from which the sample was obtained. It should also be emphasized that the single values obtained using the QC assay, do not represent the potency of the cells. Cell potency is performed using a different methodology to that used for quality control. Please see the related assays at the end of this manual.

13. HALO® QC Assay Measurement Assurance and Validation Parameters

If HALO® has been calibrated and standardized, ATP bioluminescence technology allows the User's results to be compared to the measurement assurance parameters shown in Table 4. For each control, ATP standard dose and the log-log linear regression curve fit parameters provided, the User's results must lie within the ranges provided. If this is the case, then the following are applicable:

1. The User has performed and passed the integrated proficiency test.
2. The instrument and assay readout reagents are working correctly.
3. The User can continue to process and measure samples.
4. The User can trust results of the assay.

IMPORTANT. If the User's results DO NOT comply with those in the table, DO NOT measure the samples. Perform a repeat of the controls and ATP standard curve. If the results still do not comply with those in the Table 4, contact Preferred Cell Systems for help.

Table 4: ATP Controls and Standard Curve Measurement Assurance Parameters

Expected Parameter	Observed Value	Mean \pm 15% ^(*)	Min / Max	%CV (where applicable)
0.01 μ M ATP	0.0099 μ M ATP	0.00972 - 0.0114	0.009 - 0.01	2.34%
0.03 μ M ATP	0.029 μ M ATP	0.285 - 0.0336	0.028 - 0.03	1.67%
0.05 μ M ATP	0.0497 μ M ATP	0.0486 - 0.0571	0.048 - 0.051	1.57%
0.01 μ M ATP	0.1026 μ M ATP	0.1003 - 0.118	0.099 - 0.107	1.96%
0.3 μ M ATP	0.317 μ M ATP	0.310 - 0.364	0.302 - 0.325	1.51%
0.5 μ M ATP	0.5023 μ M ATP	0.491 - 0.578	0.491 - 0.515	1.19%
1.0 μ M ATP	1.048 μ M ATP	1.024 - 1.205	0.977 - 1.117	3.7%
3.0 μ M ATP	2.722 μ M ATP	2.661 - 3.130	2.633 - 2.934	2.09%
Intercept	6.533	6.386 - 7.513	5.86 - 6.7	1.84%
Slope	0.9656	0.944 - 1.110	0.947 - 0.988	1.21%
r ² goodness of fit)	0.9993	-	0.998 - 1	0.05%
R (correlation coefficient)	1	-	0.999 - 1	0.02%

Expected Parameter	Observed Value	Mean \pm 15% ^(*)	Min / Max	%CV (where applicable)
Low control, (0.05 μ M ATP)	0.0487 μ M ATP	0.0476 - 0.0560	0.042 - 0.063	6.79%
High control 0.7 μ M ATP	0.725	0.710 - 0.836	0.655 - 0.904	5.35%
Extra high control (1.75 μ M ATP)	1.756	1.717 - 2.019	1.61 - 2.198	5.24%

(*) 15% represents the acceptable range of values for FDA Bioanalytical Method Validation Guidelines

Samples Values:

- Lowest ATP value indicating unsustainable cell proliferation for hematopoietic cells: \sim 0.04 μ M. The proliferation ability (growth) of the cells in this situation is extremely questionable.
- ATP concentration below which cells are not metabolically viable, i.e. the cells are dead: \sim 0.01 μ M. The cells should not be used in this situation.
- All sample values must lie on the ATP standard curve for accurate RLU to ATP conversion. If ATP values are greater than 3 μ M, the replicate samples should be diluted with IMDM and re-measured. Take the dilution value into account when estimating the true ATP concentration. Alternatively, repeat the culture and ATP measurement using fewer cells.

Assay Validation Parameters

HALO® exhibits the following validation parameters:

- Assay ATP linearity => 4 logs
- Assay ATP sensitivity: \sim 0.001 μ M
- Assay cell sensitivity: 20-25 cells/well (depending on cell type and purity)
- Accuracy (% correct outcomes): \sim 95%
- Sensitivity and specificity detected by Receiver Operator Characteristics (ROC) curve fit and detected as area under the curve (AUC): 0.73 - 0.752 (lowest possible value, 0.5; highest possible value, 1).
- Precision (Reliability and Reproducibility) =< 15%. At lower limit of quantification (LLOQ): 20%
- Robustness (intra- and inter-laboratory comparison): \sim 95%.
- High throughput capability (Z-Factor): >0.76 (lowest possible value, 0.5; highest possible value, 1).

14. Troubleshooting

If Calibration and Standardization Results Do Not Conform to Measurement Assurance Parameters (Section 13)

If the investigator has elected to calibrate and standardize the assay using the ATP controls and standard supplied with the kit, the results should be within the ranges provided in Section 13. If the values obtained conform to the measurement assurance parameters, the investigator can continue the assay and process and measure the samples with the assurance that the results can be trusted.

Preferred Cell Systems™

If any of the values obtained during calibration and standardization do not conform or are not within the ranges provided in Section 12, the user should repeat the calibration and standardization. Often discrepancies occur due to pipetting and/or dilution errors. Accurate and careful dilution of the ATP stock solution is important. It is also possible that if pipettes have not been professionally calibrated, errors can occur. These will also be picked up during this phase of the assay. Finally, if the ATP-ER has not been handled or stored correctly, it will decay, leading to erroneous results. Please contact Preferred Cell Systems™ to obtain new ATP-ER.

High Coefficients of Variation (%CV)

Coefficients of variation (%CV) should be $\leq 15\%$. The percent coefficient of variation is calculated as standard deviation/mean $\times 100$. High %CVs are usually an indication of incorrect dilutions or pipetting error. Although outliers can be obtained, these being observed for the more primitive stem cells than for the more mature proliferating cells, large variations between replicates should not be obtained. Please consider the following:

- Accurate reagent dispensing and mixing are of prime importance. Since the volumes dispensed are small it is imperative to use instruments that have been properly calibrated to avoid pipetting error.
- Insufficient mixing of components prior to cell plating and insufficient mixing during the addition of luminescence reagents to cultures in the 96-well plate can also lead to high CVs. Use repeater pipettes. Use calibrated or self-calibrating electronic pipettes or dispensers to add and mix the luminescence reagents.
- If the luminometer requires determining the “gain” empirically, it is possible that this parameter has not been optimally set and will result in an incorrect signal to noise ratio. Once the optimal “gain” has been set for the instrument, it should not be changed.

Low RLU Values

Performing an ATP dose response prior to sample measurement can help detect problems prior to sample measurement. If low RLU values occur, this could be due to the following reasons.

- *Reagent decay:* The ATP-ER decays with time, even when frozen. This can lead to low bioluminescence. Once thawed the reagent can be refrozen up to 11 cycles without significant loss of sensitivity. Do not use the reagent after expiry date has elapsed. As a rule of thumb, the RLU value for the lowest ATP standard should be 10 times greater than that of the background value.
- *Inadequate cell growth:* Cells did not exhibit sufficiently high viability. Measure cell viability prior to using cells. A cell viability lower than 85% should not be used. Viabilities lower than 85% can be an indication that the sample was not processed in a time-sensitive manner or that the processing procedures were not standardized and controlled. **PLEASE NOTE:** Dye exclusion viability can produce a false positive results. Even if the cell viability by dye exclusion methods is greater than 85%, this does not mean that the metabolic viability allows the cells to proliferate and grow. It is possible to have high dye exclusion viability values, but the proliferation ability of the cells will be low or non-existent, i.e. the cells might be dead.
- *Reagent deterioration:* Reagents arrived thawed, at room temperature or greater or were not stored correctly.
- *Inadequate incubator conditions:* Maintaining a correct humidified gaseous atmosphere in the incubator is essential (See Culture Plate Drying Out).
- *Carbon dioxide concentration is inadequate:* Ensure that the carbon dioxide concentration in the incubator is correct using a Fyrite gas analyzer.
- *Use low oxygen tension:* Using an oxygen concentration of 5% reduces oxygen toxicity due to free radical production and increases plating efficiency. Check that the incubator oxygen concentration is correct using a Fyrite gas analyzer.
- *Low humidity:* Plates dry out (see below) and cell growth declines.

Preferred Cell Systems™

- *Contamination:* Cells cultured in 96-well plates cannot be view under a microscope. If contamination occurs it will usually be seen by the difference in color of the cultures, if the medium contains an indicator, e.g. phenol red. Contaminated cultures will usually be bright yellow in color and probably cloudy in appearance. Cell cultures that demonstrate high proliferation will usually appear orange to light orange, but will not be cloudy. If only “spot” contamination occurs, this is usually due to pipette or repeater tips coming in contact with materials other than the reagents. Contamination will usually lead to outlier RLU values.

Luminescence Reagent Mixing.

The luminescence reagent has to be added and thoroughly mixed with the culture components. The ATP-ER lyses the cells and releases intracellular ATP. If mixing is not adequate, only a proportion of the cells will be lysed and the RLU values will be low. Conversely, too much mixing can lead to ATP degradation and low luminescence readings.

Culture Plates Drying Out

- Due to the relatively small culture volume (0.1mL), drying out of the culture wells, particularly around the outside of the plate can be a problem. These are called “edge effects”. An incubator with insufficient humidity will cause this problem. To ensure that this does not occur, the incubator water reservoir should be full and the humidity in the chamber checked using a hygrometer.
- If drying out continues, use of a humidity chamber is recommended. Please refer to Section 9 (v) for instructions on how to build a humidity chamber.

Other Hematopoietic Cellular Therapy Assay Products from Preferred Cell Systems™

HALO® PCA is a 96-well ATP bioluminescence progenitor cell assay to replace the CFU assay for cell processing laboratories.

HemoFLUOR™ PCA is a 96-well, fluorescence readout version of HALO® QC.

HemoLIGHT™ PCA is a 96-well, absorbance readout version of HALO® QC.

STEMpredict™ is a 3-day, fully standardized, ATP bioluminescence, stem cell quality assay designed primarily for cord blood banks to triage high from low quality cord blood units prior to cryopreservation. STEMpredict™ is the only assay for hematopoietic cellular therapy products designed for both high-throughput 96- and 384-well plate formats.

HALO® RS is a 7 day, 96-well, standardized ATP bioluminescence assay to establish cord blood, bone marrow or peripheral blood reference standards for HALO® SC-IPS assays.

HALO® SC-IPS is a 7-day, standardized, 96-well, ATP bioluminescence assay to measure the identity, purity and strength (potency) of 2 primitive stem cell populations in cord blood, bone marrow or peripheral blood samples, prior to use in patients.

HALO® TE, HemoFLUOR™ TE and HemoLIGHT™ TE are time to engraftment assays.

HALO® PMT, HemoFLUOR™ PMT and HemoLIGHT™ PMT are “global” 4-, 5- or 7-population assays to monitor patient reconstitution after transplantation.

Ordering Information

Toll free: 1-888-436-6869

Tel: (719) 264-6251

Fax: (719) 264-6253

Email: info@preferred-cell-systems.com

Order online at preferred-cell-systems.com

Technical Support

Tel: (719) 264-6251

Email: info@preferred-cell-systems.com

Preferred Cell Systems™

1485 Garden of the Gods Road

Suite 152

Colorado Springs, CO 80907

U.S.A.


Website: www.preferred-cell-systems.com

HALO® is a registered trademark of HemoGenix®, Inc
HALO® was designed and developed by HemoGenix®, Inc.
Assays-in-a-Box™ and Preferred Cell Systems™ are trademarks of Preferred Cell Systems™

Preferred Cell Systems™

Calibration and Standardization Protocol of an ATP Bioluminescence Assay

PROTOCOL 1: ATP Standard Curve from 0.01 μ M to 1 μ M For Samples with Known or Expected Normal Cell Proliferation



Calibration and Standardization Protocol of an ATP Bioluminescence Assay

PROTOCOL 2: ATP Standard Curve from 0.03 μ M - 3 μ M For Samples with Known or Expected High Cell Proliferation

