

HemoFLUORTM PCA

— By Preferred Cell SystemsTM —


Fluorescence Progenitor Cell Assay Methylcellulose-Free CFU Replacement for Cell Processing Laboratories

Uses and Benefits

- Replaces all methylcellulose based reagents. No colony counting necessary.
- Fluorescence readout combined with SECTM Technology for the highest precision.
- Short 4-6 culture depending on species. Results in 30 min. or less.
- For all transplantation tissue.
- Highly reproducible results.
- Complete assay kits. Reagents and 96-well plate included for up to 24 samples/plate.
- Smaller sample and reagent volume. Faster setup.
- Optional standardization and validation available.
- Simple to both learn and use.
- Similar growth factor cocktails to MethoCultTM and other methylcellulose reagents, but far superior cell growth. Direct correlation with CFU.
- Detect primitive stem and progenitor cells easily and efficiently.
- Saves time and labor allowing personnel to perform other important activities.

Assay Principle

Incorporates the measurement of non-fluorescent Resazurin (AlimaBlueTM) to fluorescent resorufin in an oxidation/reduction reaction that occurs only in live cells. The amount of fluorescence is dependent on the metabolic status of the cells and, therefore, is a measure of metabolic viability and cell proliferation or inhibition. After adding the reagent, the signal develops within 4 hours in a fluorescence plate reader using a 560nm excitation filter and a 590nm emission filter.



Assays-in-a-BoxTM
by
Preferred Cell SystemsTM

www.preferred-cell-systems.com

HemoFLUORTMPCA: A Replacement for Methylcellulose CFU Assays

HemoFLUOR TM PCA Number	HemoFLUOR TM PCA Cell Population	Equivalent CFU Cell Population	Equivalent MethoCult TM Reagent	Growth Factor Cocktail
PCA1	SC-GEM 3	CFC-GEM 3	H4434 'Classic'	EPO, GM-CSF, IL-3, SCF
PCA2	SC-GEM 2	CFC-GEM 2	H4034 'Optimum'	EPO, GM-CSF, G-CSF, IL-3, SCF
PCA3	P-GM 1	GM-CFC 1	H4534 'Classic'	GM-CSF, IL-3, SCF
PCA4	P-GM 2	GM-CFC 2	H4035 'Optimum'	GM-CSF, G-CSF, IL-3, SCF
PCA5	SC-GEMM 3	CFC-GEMM 3	H4435 'Enriched'	EPO, GM-CSF, G-CSF, IL-3, IL-6, SCF TPO(#)
PCA6	SC-GEMM	CFC-GEMM	Methocult 'Express'	Recombinant cytokines EPO

(#) denotes that the Preferred Cell SystemsTM formulation of this product includes TPO that stimulates the production of megakaryopoiesis. MethoCultTM and other competitor formulations of this specific product do not include TPO, do not stimulate the production of megakaryocytes and therefore do not detect the CFC-GEMM stem cell population.