

Certificate of Analysis

Product Name: **ATP Enumeration Reagent Lyophilized Substrate**
Lot Number: **391761**


Storage Conditions: **-20°C**
Expiration Date: **10/20/22**

FOR RESEARCH USE ONLY
CAUTION: Not intended for human or animal therapeutic use.

Quality Control

Test	Specification	Results
Signal Stability	Signal ½ life > 5hr	Pass
Sensitivity	$RLU_{0.01\mu MATP} > RLU_{Bkg} + 2(S.D._{0.01\mu MATP} + S.D._{Bkg})$	Pass
Batch to Batch Consistency	Comparison to a Dynex reference plate	Pass

Manufactured and tested for Preferred Cell Systems™ by Supplier


Paul Bork, QA Manager
Date: *3/12/2020*