

## Certificate of Analysis

Product Name: **ATP Standard, 10µM**

Lot Number: **ATP 250618**

Storage Conditions: **-20°C**

Expiration Date: **06/25/23**

**FOR RESEARCH USE ONLY**

**CAUTION: Not intended for human or animal therapeutic use.**

### **Quality Control**

Test	Specification	Results
ATP Concentration (10µM) by UV-VIS Spectrophotometry	9.0 – 11.0 µM	10.56 µM
Diluted ATP Concentration (1µM) by Gold Standard ATP Curve	0.900 – 1.100 µM	1.021 µM

Manufactured and tested at Preferred Cell Systems™



Paul Bork, QA Manager

Date: *7/23/19*